

EasyShift™



Patient Repositioning System

Turn and Boost Patients at the Push of a Button up to 500 lbs.

EasyShift™ Overlay is installed over your current mattress surface.

Air Pillows inflate on left, right, and top sides of the patient.

Top Air Pillow shifts the patient up in bed.

Left and Right Air Pillows laterally shift the patient during the turn to avoid hitting the siderail.

- Industry-leading 90° rotation
- Air Pillows help hold patient in 30° offloading position

Cleaned just like a regular mattress cover. No special sheets or disposables.



Patient Benefits

- Frequent turning is the best way to reduce pressure injuries.
- Air Pillows inflate gently and gradually, reducing rapid shear forces that can cause skin abrasions and lacerations.
- Keeping patient elevated at 30° improves respiratory function.

Caregiver Benefits

- Streamlines workflow for routine care activities. Patients of size can be repositioned by 1 staff member, saving time for other caregivers.
- EasyShift™ virtually eliminates the #1 and #2 causes of lifting strain in acute care.